

ZNF280D Rabbit Polyclonal Antibody



CAB16538

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

17kDa/24kDa/48kDa/65kDa/
107kDa/109kDa

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

54816

Uniprot

Q6N043

Synonyms:

ZNF280D; SUHW4; ZNF634

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 700-979 of human ZNF280D (NP_060131.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
