

ANO1 Rabbit Polyclonal Antibody



CAB16540

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

74kDa/97kDa/114kDa

Applications:

IF

Reactivity:

Human

Antibody Information

Recommended dilutions:

IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

55107

Uniprot

Q5XXA6

Synonyms:

ANO1; DOG1; ORAOV2; TAOS2; TMEM16A

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 40-210 of human ANO1 (NP_060513.5).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
