

RHBG Rabbit Polyclonal Antibody



CAB16558

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

32kDa/39kDa/42kDa/43kDa/
47kDa

Applications:

WB

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes one of two non-erythroid members of the Rhesus (Rh) protein family. Non-erythroid Rh protein family members are mainly expressed in the kidney and belong to the methylammonium-ammonium permease/ammonia transporters superfamily. All Rh family proteins are predicted to be transmembrane proteins with 12 membrane spanning domains and intracytoplasmic N- and C-termini. Alternative splicing of this gene results in multiple transcript variants.

Immunogen information

Gene ID:

57127

Uniprot

Q9H310

Synonyms:

RHBG; SLC42A2

Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RHBG (NP_065140.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
