

KIF17 Rabbit Polyclonal Antibody



CAB16562

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

170kDa

Calculated MW:

114kDa/115kDa

Applications:

WB IHC IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

57576

Uniprot

Q9P2E2

Synonyms:

KIF17; KIF17B; KIF3X; KLP-2; OSM-3


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 360-660 of human KIF17 (NP_001116291).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images




Western blot analysis of extracts of various cell lines, using KIF17 antibody (CAB16562) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 30s.



Immunohistochemistry of paraffin-embedded Human lung cancer using KIF17 Rabbit pAb (CAB16562) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Rat testis using KIF17 Rabbit pAb (CAB16562) at dilution of 1:100 (40x lens).



Immunofluorescence analysis of rat brain using KIF17 antibody (CAB16562) at dilution of 1:100. Blue: DAPI for nuclear staining.