

Product Information

Product SKU:	CAB16570	Gene ID:	79567	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	140kDa	Conjugate:	Unconjugated
Calculated MW:	132kDa	Isotype:	IgG

Immunogen Information

Background: Enables 14-3-3 protein binding activity. Involved in several processes, including establishment of Golgi localization; negative regulation of Rho guanyl-nucleotide exchange factor activity; and negative regulation of Rho protein signal transduction. Located in Golgi apparatus; cell leading edge; and membrane.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: FAM65A

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1004-1223 of human FAM65A (NP_001180451.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.