

ZGPAT Rabbit Polyclonal Antibody



CAB16583

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

57kDa

Calculated MW:

20kDa/54kDa/55kDa/57kDa

Applications:

WB

Reactivity:

Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

84619

Uniprot

Q8N5A5

Synonyms:

ZGPAT; GPATC6; GPATCH6; KIAA1847; ZC3H9; ZC3HDC9; ZIP


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 212-511 of human ZGPAT (NP_852150.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of rat lung, using ZGPAT antibody (CAB16583) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 3s.