

Product Information

Product SKU:	CAB16586	Gene ID:	84967	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	14kDa	Conjugate:	Unconjugated
Calculated MW:	14kDa	Isotype:	IgG

Immunogen Information

Background: Enables U7 snRNA binding activity. Involved in positive regulation of G1/S transition of mitotic cell cycle. Located in Cajal body. Part of U7 snRNP.

Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	MST074; MSTP074; LSM10
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-123 of human LSM10 (NP_116270.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.