

Product Information

Product SKU:	CAB16587	Gene ID:	89839	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	30kDa	Conjugate:	Unconjugated
Calculated MW:	30kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable GTPase activator activity. Involved in cerebral cortex development and negative regulation of mitochondrial membrane permeability. Acts upstream of with a positive effect on glutamine catabolic process. Located in mitochondrial matrix.

Recommended Dilution:	WB,1:100 - 1:500
Synonyms:	B'-T; FAM7B1; GAP (1-8); ARHGAP11B
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 218-267 of human ARHGAP11B (NP_001034930.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.