

# SYNPR Rabbit Polyclonal Antibody



CAB16595

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

29kDa/31kDa

### Applications:

WB

### Reactivity:

Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

132204

#### Uniprot

Q8TBG9

#### Synonyms:

SYNPR; SPO

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 196-265 of human SYNPR (NP\_653243.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---