

# GPR142 Rabbit Polyclonal Antibody



**CAB16617**

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

51kDa

### Applications:

WB

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

GPR142 is a member of the rhodopsin family of G protein-coupled receptors (GPRs).

## Immunogen information

### Gene ID:

350383

### Uniprot

Q7Z601

### Synonyms:

GPR142; GPRg1b; PGR2

### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human GPR142 (NP\_861455.1).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---