

# HRNR Rabbit Polyclonal Antibody



CAB16620

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

### Calculated MW:

282kDa

### Applications:

IF

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

IF 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

388697

#### Uniprot

Q86YZ3

#### Synonyms:

HRNR; FLG3; S100A16; S100a18; hornerin

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human HRNR (NP\_001009931.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunofluorescence analysis of HeLa cells using HRNR antibody (CAB16620) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.