

CAB16635

Product Information

Product SKU:	CAB16635	Gene ID:	347252	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	35kDa	Conjugate:	Unconjugated
Calculated MW:	29kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable insulin-like growth factor binding activity. Involved in cellular response to tumor cell. Located in extracellular space.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	IGFBPRP4; IGFBP-RP4; bA113O24.1; IGFBPL1
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human IGFBPL1 (NP_001007564.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.