

Product Information

Product SKU:	CAB16642	Gene ID:	6568	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Rat

Additional Information

Observed MW:	51kDa	Conjugate:	Unconjugated
Calculated MW:	51kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable sialic acid transmembrane transporter activity. Involved in urate metabolic process and urate transport. Located in apical plasma membrane.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	NPT1; NPT-1; NAPI-1; SLC17A1
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SLC17A1 (NP_005065.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.