

ZDHHC14 Rabbit Polyclonal Antibody



CAB16644

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

53kDa

Calculated MW:

51kDa/53kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

224454

Uniprot

Q8BQQ1

Synonyms:

ZDHHC14; NEW1CP


Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse ZDHHC14 (NP_666185.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - ZDHHC14 Rabbit pAb (CAB16644)