

Product Information

Product SKU:	CAB16653	Gene ID:	55273	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	14kDa	Conjugate:	Unconjugated
Calculated MW:	14kDa	Isotype:	IgG

Immunogen Information

Background: Involved in BMP signaling pathway. Located in several cellular components, including endoplasmic reticulum; perikaryon; and perinuclear region of cytoplasm.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	TMEM100
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 50 to the C-terminus of human TMEM100 (NP_001093110.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.