

CALN1 Rabbit Polyclonal Antibody



CAB16684

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

24kDa/29kDa

Applications:

WB

Reactivity:

Human

Protein Background

This gene encodes a protein with high similarity to the calcium-binding proteins of the calmodulin family. The encoded protein contains two EF-hand domains and potential calcium-binding sites. Alternative splicing results in multiple transcript variants.

Immunogen information

Gene ID:

83698

Uniprot

Q9BXU9

Synonyms:

CALN1; CABP8

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 100-180 of human CALN1 (NP_001017440.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
