

Product Information

Product SKU:	CAB16726	Gene ID:	5875	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	55kDa	Conjugate:	-
Calculated MW:	65kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable small GTPase binding activity. Predicted to contribute to Rab geranylgeranyltransferase activity. Predicted to be involved in protein geranylgeranylation. Located in nucleoplasm.

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Synonyms:	PTAR3; RABGGTA
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 240-360 of human RABGGTA (NP_878256.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.