

Product Information

Product SKU:	CAB16741	Gene ID:	270328	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	54kDa	Conjugate:	Unconjugated
Calculated MW:	54kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable phosphatidylinositol-4,5-bisphosphate binding activity; phosphatidylinositol-4-phosphate binding activity; and phosphatidylserine binding activity. Predicted to be involved in defense response to bacterium and pyroptosis. Predicted to act upstream of or within programmed cell death. Predicted to be located in cytoplasm. Is expressed in gut and metanephros. Orthologous to human GSDMC (gasdermin C).

Recommended Dilution:	WB,1:1000 - 1:5000
Synonyms:	9930109F21Rik; Gsdmc3
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 2-80 of mouse Gsdmc3 (NP_899017.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.