

Gsdmc3 Rabbit Polyclonal Antibody



CAB16741

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

54KDa

Calculated MW:

54kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

270328

Uniprot

Q8CB12

Synonyms:

GSDMC; MLZE; 9930109F21Rik; gasdermin-C; Gsdmc3

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 2-80 of mouse Gsdmc3 (NP_899017.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - Gsdmc3 Rabbit pAb (CAB16741)

