

# APOF Rabbit Polyclonal Antibody



**CAB16753**

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

35kDa

### Applications:

WB

### Reactivity:

Rat

## Protein Background

The product of this gene is one of the minor apolipoproteins found in plasma. This protein forms complexes with lipoproteins and may be involved in transport and/or esterification of cholesterol.

## Immunogen information

### Gene ID:

319

### Uniprot

Q13790

### Synonyms:

APOF; Apo-F; LTIP

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:1000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Immunogen:

A synthetic Peptide of human APOF

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

## Product Images

---