

CAB16762

Product Information

Product SKU:	CAB16762	Gene ID:	388	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	22kDa	Conjugate:	Unconjugated
Calculated MW:	22kDa	Isotype:	IgG

Immunogen Information

Background:	Predicted to enable GTP binding activity; GTPase activity; and protein kinase binding activity. Involved in several processes, including cellular response to hydrogen peroxide; cellular response to ionizing radiation; and regulation of cell migration. Located in cleavage furrow and endosome membrane. Biomarker of breast cancer.
Recommended Dilution:	WB,1:500 - 1:2000 IP,0.5µg-4µg antibody for 200µg-400µg extracts of whole cells
Synonyms:	ARH6; ARHB; RHOH6; MST081; MSTP081; RhoB
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 100-196 of human RhoB (NP_004031.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.