

BBS1 Rabbit Polyclonal Antibody



CAB16773

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

49kDa/65kDa/69kDa

Applications:

IHC

Reactivity:

Human, Mouse

Protein Background

Mutations in this gene have been observed in patients with the major form (type 1) of Bardet-Biedl syndrome. The encoded protein may play a role in eye, limb, cardiac and reproductive system development.

Immunogen information

Gene ID:

582

Uniprot

Q8NFJ9

Synonyms:

BBS1; BBS2L2

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 240-370 of human BBS1 (NP_078925.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
