

# C5AR1 Rabbit Polyclonal Antibody



**CAB16784**

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

39kDa

### Applications:

IHC

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

IHC 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

728

#### Uniprot

P21730

#### Synonyms:

C5AR1; C5A; C5AR; C5R1; CD88

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 250-350 of human C5AR1 (NP\_001727.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---