

CD52 Rabbit Polyclonal Antibody



CAB16815

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

29kDa

Calculated MW:

Applications:

WB IHC IF

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:100 IF 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

1043

Uniprot

P31358

Synonyms:

CD52; CDW52; EDDM5

Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 1-61 of human CD52 (NP_001794.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - CD52 Rabbit pAb (CAB16815)

