

Product Information

Product SKU:	CAB16815	Gene ID:	1043	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	29kDa	Conjugate:	-
Calculated MW:	7kDa	Isotype:	IgG

Immunogen Information

**Background:** Involved in positive regulation of cytosolic calcium ion concentration. Predicted to be located in extracellular region and plasma membrane. Predicted to be intrinsic component of plasma membrane. Predicted to be active in sperm midpiece.

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:100
Synonyms:	HE5; CDW52; EDDM5; CD52
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 1-61 of human CD52 (NP_001794.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.