

---

## Product Information

<b>Product SKU:</b>	CAB16845	<b>Gene ID:</b>	2114	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	Refer to figures	<b>Conjugate:</b>	-
<b>Calculated MW:</b>	53kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** This gene encodes a transcription factor which regulates genes involved in development and apoptosis. The encoded protein is also a protooncogene and shown to be involved in regulation of telomerase. A pseudogene of this gene is located on the X chromosome. Alternative splicing results in multiple transcript variants.

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** ETS2IT1; ETS2

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 300-400 of human ETS2 (NP\_005230.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.