

FXN Rabbit Polyclonal Antibody



CAB16853

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

26kDa

Calculated MW:

19kDa/21kDa/23kDa

Applications:

WB IHC IF

Reactivity:

Human, Mouse, Rat

Protein Background

This nuclear gene encodes a mitochondrial protein which belongs to the FRATAXIN family. The protein functions in regulating mitochondrial iron transport and respiration. The expansion of intronic trinucleotide repeat GAA from 8-33 repeats to >90 repeats results in Friedreich ataxia. Alternative splicing results in multiple transcript variants.

Immunogen information

Gene ID:

2395

Uniprot

Q16595

Synonyms:

FXN; CyaY; FA; FARR; FRDA; X25; frataxin

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 42-210 of human FXN (NP_000135.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - FXN Rabbit pAb (CAB16853)