

LRRC32 Rabbit Polyclonal Antibody



CAB16858

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a type I membrane protein which contains 20 leucine-rich repeats. Alterations in the chromosomal region 11q13-11q14 are involved in several pathologies.

Immunogen information

Gene ID:

2615

Uniprot

Q14392

Synonyms:

LRRC32; D11S833E; GARP

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 520-620 of human LRRC32 (NP_005503.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
