
Product Information

Product SKU:	CAB16888	Gene ID:	23362	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	116kDa	Conjugate:	-
Calculated MW:	116kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable guanyl-nucleotide exchange factor activity and phospholipid binding activity. Predicted to be involved in regulation of ARF protein signal transduction and regulation of catalytic activity. Predicted to be located in membrane. Predicted to be active in ruffle membrane.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: EFA6D; EFA6R; HCA67; PSD3

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 878-1047 of human PSD3 (NP_056125.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.