

Product Information

Product SKU:	CAB16925	Gene ID:	4707	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	11kDa	Conjugate:	-
Calculated MW:	7kDa	Isotype:	IgG

Immunogen Information

Background: Involved in mitochondrial respiratory chain complex I assembly. Located in mitochondrion and nuclear speck. Part of mitochondrial respiratory chain complex I.

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:100
Synonyms:	MNLL; CI-MNLL; CI-SGDH; NDUFB1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human NDUFB1 (NP_004536.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.