

NPAT Rabbit Polyclonal Antibody



CAB16931

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

4863

Uniprot

Q14207

Synonyms:

NPAT; E14; E14/NPAT; p220

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1128-1427 of human NPAT (NP_002510.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
