

Product Information

Product SKU:	CAB16934	Gene ID:	4887	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	42kDa	Conjugate:	-
Calculated MW:	43kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable calcium channel regulator activity and neuropeptide Y receptor activity. Involved in cardiac left ventricle morphogenesis and outflow tract morphogenesis. Located in cilium. Implicated in Huntington's disease; morbid obesity; and obesity. Biomarker of peripheral artery disease and temporal lobe epilepsy.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	NPY2-R; NPY2R
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 300 to the C-terminus of human NPY2R (NP_000901.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.