

CAB16985

Product Information

Product SKU:	CAB16985	Gene ID:	7057	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	Refer to figures	Conjugate:	-
Calculated MW:	129kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is a subunit of a disulfide-linked homotrimeric protein. This protein is an adhesive glycoprotein that mediates cell-to-cell and cell-to-matrix interactions. This protein can bind to fibrinogen, fibronectin, laminin, type V collagen and integrins alpha-V/beta-1. This protein has been shown to play roles in platelet aggregation, angiogenesis, and tumorigenesis.

Recommended Dilution:	IHC-P, 1:50 - 1:200
Synonyms:	TSP; THBS; TSP1; TSP-1; THBS-1; THBS1
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human THBS1 (NP_003237.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.