

# TSPY1 Rabbit Polyclonal Antibody



CAB16992

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

### Applications:

WB

### Reactivity:

Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

The protein encoded by this gene is found only in testicular tissue and may be involved in spermatogenesis. Many functional paralogs and pseudogenes of this gene are present in a cluster in humans, but only a single, nonfunctional orthologous gene is found in mouse. Alternative splicing results in multiple transcript variants.

## Immunogen information

### Gene ID:

7258

### Uniprot

Q01534

### Synonyms:

TSPY1; CT78; DYS14; TSPY; pJA923

### Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TSPY1 (NP\_001071165.2).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---