
Product Information

Product SKU:	CAB17042	Gene ID:	9024	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	78-90kDa	Conjugate:	-
Calculated MW:	82kDa	Isotype:	IgG

Immunogen Information

Background:	Enables several functions, including ATP binding activity; ATPase binding activity; and protein kinase activity. Involved in several processes, including cellular protein metabolic process; intrinsic apoptotic signaling pathway in response to endoplasmic reticulum stress; and regulation of insulin secretion involved in cellular response to glucose stimulus. Located in centrosome and endoplasmic reticulum.
Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	SAD1; SADA; SAD-A; STK29; PEN11B; C11orf7; BRSK2
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human BRSK2 (NP_001243556.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.