

# BRSK2 Rabbit Polyclonal Antibody



CAB17042

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

68kDa/74kDa/75kDa/77kDa/  
81kDa/84kDa

### Applications:

WB

### Reactivity:

Human, Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

9024

#### Uniprot

Q8IWQ3

#### Synonyms:

BRSK2; C11orf7; PEN11B; SAD1; SADA; STK29

#### Immunogen:

A synthetic peptide of human BRSK2.

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---