

TAOK2 Rabbit Polyclonal Antibody



CAB17048

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

138kDa

Calculated MW:

118kDa/119kDa/126kDa/138kDa

Applications:

WB

Reactivity:

Human

Protein Background

This gene encodes a serine/threonine protein kinase that is involved in many different processes, including, cell signaling, microtubule organization and stability, and apoptosis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

Immunogen information

Gene ID:

9344

Uniprot

Q9UL54

Synonyms:

TAOK2; MAP3K17; PSK; PSK1; PSK1-BETA; TAO1; TAO2

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 610-730 of human TAOK2 (NP_057235.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - TAOK2 Rabbit pAb (CAB17048)