

# RAB33A Rabbit Polyclonal Antibody



CAB17049

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

27kDa

### Calculated MW:

### Applications:

WB

### Reactivity:

Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

The protein encoded by this gene belongs to the small GTPase superfamily, Rab family. It is GTP-binding protein and may be involved in vesicle transport.

## Immunogen information

### Gene ID:

9363

### Uniprot

Q14088

### Synonyms:

RAB33A; RabS10

### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 118-237 of human RAB33A (NP\_004785.1).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

Western blot - RAB33A Rabbit pAb (CAB17049)

