
Product Information

Product SKU:	CAB17049	Gene ID:	9363	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	27kDa	Conjugate:	-
Calculated MW:	27kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene belongs to the small GTPase superfamily, Rab family. It is GTP-binding protein and may be involved in vesicle transport.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: RabS10; RAB33A

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 118-237 of human RAB33A (NP_004785.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.