

LRRC14 Rabbit Polyclonal Antibody



CAB17057

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

54kDa

Calculated MW:

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a leucine-rich repeat-containing protein. Alternate splicing results in multiple transcript variants.

Immunogen information

Gene ID:

9684

Uniprot

Q15048

Synonyms:

LRRC14; LRRC14A


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human LRRC14 (NP_055480.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - LRRC14 Rabbit pAb (CAB17057)