

VGLL4 Rabbit Polyclonal Antibody



CAB17058

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

31kDa

Calculated MW:

30kDa

Applications:

WB

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

9686

Uniprot

Q14135

Synonyms:

VGLL4; VGL-4

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 96-296 of human VGLL4 (NP_001121691.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - VGLL4 Rabbit pAb (CAB17058)

