

RAB35 Rabbit Polyclonal Antibody



CAB17091

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

17kDa/23kDa

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

11021

Uniprot

Q15286

Synonyms:

RAB35; H-ray; RAB1C; RAY

Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 50-150 of human RAB35 (NP_006852.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
