

# NID2 Rabbit Polyclonal Antibody



CAB17097

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

### Applications:

WB

### Reactivity:

Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

This gene encodes a member of the nidogen family of basement membrane proteins. This protein is a cell-adhesion protein that binds collagens I and IV and laminin and may be involved in maintaining the structure of the basement membrane.

## Immunogen information

### Gene ID:

22795

### Uniprot

Q14112

### Synonyms:

NID2; NID-2; nidogen-2

### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 900-1200 of human NID2 (NP\_031387.3).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---