

NID2 Rabbit Polyclonal Antibody



CAB17098

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

This gene encodes a member of the nidogen family of basement membrane proteins. This protein is a cell-adhesion protein that binds collagens I and IV and laminin and may be involved in maintaining the structure of the basement membrane.

Immunogen information

Gene ID:

22795

Uniprot

Q14112

Synonyms:

NID2; NID-2; nidogen-2

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 900-1200 of human NID2 (NP_031387.3).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
