

CBX7 Rabbit Polyclonal Antibody



CAB17115

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

IHC

Reactivity:

Human, Mouse, Rat

Protein Background

This gene encodes a protein that contains the CHROMO (CHRomatin Organization MODifier) domain. The encoded protein is a component of the Polycomb repressive complex 1 (PRC1), and is thought to control the lifespan of several normal human cells.

Immunogen information

Gene ID:

23492

Uniprot

O95931

Synonyms:

CBX7

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human CBX7 (NP_783640.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
