

# PPP1R15A Rabbit Polyclonal Antibody



CAB17117

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

49kDa/50kDa 56kDa/73kDa

### Applications:

WB

### Reactivity:

Human, Mouse, Rat

## Protein Background

This gene is a member of a group of genes whose transcript levels are increased following stressful growth arrest conditions and treatment with DNA-damaging agents. The induction of this gene by ionizing radiation occurs in certain cell lines regardless of p53 status, and its protein response is correlated with apoptosis following ionizing radiation.

## Immunogen information

### Gene ID:

23645

### Uniprot

O75807

### Synonyms:

PPP1R15A; GADD34

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 375-674 of human PPP1R15A (NP\_055145.3).

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---