

ARFIP2 Rabbit Polyclonal Antibody



CAB17118

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

29kDa/33kDa/37kDa

Applications:

IHC

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

23647

Uniprot

P53365

Synonyms:

ARFIP2; POR1

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human ARFIP2 (NP_036534.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
