

Product Information

| | | | | | |
|---------------------|----------|----------------------|--------|--------------------|-----------------|
| Product SKU: | CAB17123 | Gene ID: | 25820 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Mouse,Rat |

Additional Information

| | | | |
|-----------------------|-------|-------------------|-----|
| Observed MW: | 64kDa | Conjugate: | - |
| Calculated MW: | 64kDa | Isotype: | IgG |

Immunogen Information

Background: Enables enzyme binding activity; ubiquitin-protein transferase activity; and zinc ion binding activity. Involved in protein ubiquitination. Located in Lewy body; cytoplasm; and nuclear body. Colocalizes with cullin-RING ubiquitin ligase complex.

| | |
|------------------------------|--|
| Recommended Dilution: | WB,1:500 - 1:1000 |
| Synonyms: | ARI; HARI; HHARI; UBCH7BP; ARIH1 |
| Purification Method: | Affinity purification |
| Immunogen: | A synthetic peptide corresponding to a sequence within amino acids 400 to the C-terminus of human ARIH1 (NP_005735.2). |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3. |