

# ARIH1 Rabbit Polyclonal Antibody



CAB17123

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

### Applications:

WB

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

25820

#### Uniprot

Q9Y4X5

#### Synonyms:

ARIH1; ARI; HARI; HHARI; UBCH7BP

#### Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 400 to the C-terminus of human ARIH1 (NP\_005735.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---