

# GAPVD1 Rabbit Polyclonal Antibody



CAB17127

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

165kDa

### Calculated MW:

165kDa

### Applications:

WB IF

### Reactivity:

Mouse

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:100

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

26130

#### Uniprot

Q14C86

#### Synonyms:

GAPVD1; GAPEX5; GAPex-5; RAP6

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human GAPVD1 (NP\_056450.2).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot - GAPVD1 Polyclonal Antibody (CAB17127)